

## Remote Communication

1.36

©Colorimetry Research, Inc. 2018

# Contents

<b>1</b>	<b>Introduction</b>	<b>3</b>
<b>2</b>	<b>Remote Communication</b>	<b>5</b>
2.1	Installing Drivers . . . . .	5
2.1.1	Windows . . . . .	5
2.1.2	OSX . . . . .	6
2.2	Remote Communication . . . . .	6
2.2.1	Response . . . . .	7
2.2.2	Coding . . . . .	8
<b>3</b>	<b>Getting Started</b>	<b>9</b>
3.1	Getting started with the Terminal Software . . . . .	9
3.1.1	Windows® . . . . .	9
3.1.1.1	Hyper Terminal . . . . .	9
3.1.1.2	Windows®: PuTTY . . . . .	9
3.1.2	OS X . . . . .	9
3.1.2.1	screen . . . . .	9
3.2	Getting started with Remote messaging . . . . .	9
<b>4</b>	<b>Remote Commands</b>	<b>11</b>
4.1	Echo Command . . . . .	11
4.2	Setup Commands . . . . .	11
4.2.1	Measurement Setup Commands . . . . .	11
4.2.1.1	SM Accessory . . . . .	12
4.2.1.2	SM Filter1 . . . . .	12
4.2.1.3	SM Filter2 . . . . .	13
4.2.1.4	SM Filter3 . . . . .	13
4.2.1.5	SM Aperture . . . . .	14
4.2.1.6	SM Mode . . . . .	14
4.2.1.7	SM ExposureMode . . . . .	14
4.2.1.8	SM Exposure . . . . .	15
4.2.1.9	SM MaxAutoExposure . . . . .	15

4.2.1.10	SM RangeMode	16
4.2.1.11	SM Range	16
4.2.1.12	SM SyncMode	16
4.2.1.13	SM SyncFreq	17
4.2.1.14	SM ExposureX	17
4.2.1.15	SM MatrixMode	17
4.2.1.16	SM UserCalibMode	18
4.2.1.17	SM Matrix	18
4.2.1.18	SM Match	19
4.2.1.19	SM Speed	19
4.2.1.20	SM SamplingRate	20
4.2.1.21	SM MaxFreqFlickerSearch	20
4.2.1.22	SM CMF	20
4.2.1.23	SM Reset	21
4.2.2	Hardware Setup Commands	21
4.2.3	Computation Setup Commands	21
4.2.3.1	SC CMF	21
4.3	Read Commands	22
4.3.1	Read Measurement Commands	22
4.3.1.1	RM ID	23
4.3.1.2	RM Model	24
4.3.1.3	RM Time	24
4.3.1.4	RM Accessory	24
4.3.1.5	RM Filter	24
4.3.1.6	RM Aperture	25
4.3.1.7	RM Mode	25
4.3.1.8	RM ExposureMode	25
4.3.1.9	RM Exposure	26
4.3.1.10	RM MaxAutoExposure	26
4.3.1.11	RM RangeMode	26
4.3.1.12	RM Range	27
4.3.1.13	RM SyncMode	27
4.3.1.14	RM SyncFreq	27
4.3.1.15	RM ExposureX	28
4.3.1.16	RM MatrixMode	28
4.3.1.17	RM UserCalibMode	28
4.3.1.18	RM Matrix	29
4.3.1.19	RM Match	29
4.3.1.20	RM Speed	29
4.3.1.21	RM X	30

---

4.3.1.22	RM X10	30
4.3.1.23	RM Y	31
4.3.1.24	RM Y10	31
4.3.1.25	RM Z	31
4.3.1.26	RM Z10	31
4.3.1.27	RM XYZ	32
4.3.1.28	RM XYZ10	32
4.3.1.29	RM xy	32
4.3.1.30	RM xy10	33
4.3.1.31	RM uv	33
4.3.1.32	RM upvp	33
4.3.1.33	RM CCT	34
4.3.1.34	RM Warnings	34
4.3.1.35	RM Yv	34
4.3.1.36	RM Radiometric	34
4.3.1.37	RM Spectrum	35
4.3.1.38	RM Temporal	35
4.3.1.39	RM TemporalY	36
4.3.1.40	RM SamplingRate	36
4.3.1.41	RM CMF	37
4.3.2	Read Configuration Commands	37
4.3.2.1	RC ID	38
4.3.2.2	RC Model	38
4.3.2.3	RC InstrumentType	38
4.3.2.4	RC Accessory	39
4.3.2.5	RC Filter	39
4.3.2.6	RC Aperture	40
4.3.2.7	RC Mode	40
4.3.2.8	RC ExposureMode	41
4.3.2.9	RC RangeMode	41
4.3.2.10	RC Range	41
4.3.2.11	RC SyncMode	42
4.3.2.12	RC Firmware	43
4.3.2.13	RC MatrixMode	43
4.3.2.14	RC UserCalibMode	44
4.3.2.15	RC Matrix	44
4.3.2.16	RC Match	44
4.3.2.17	RC MatrixCalibration	45
4.3.2.18	RC MatrixCalib	45
4.3.2.19	RC MatchCalib	45

---

4.3.2.20	RC MinExposure	46
4.3.2.21	RC MaxExposure	46
4.3.2.22	RC MinSyncFreq	46
4.3.2.23	RC MaxSyncFreq	47
4.3.2.24	RC MinExposureX	47
4.3.2.25	RC MaxExposureX	48
4.3.2.26	RC Speed	48
4.3.2.27	RC MinSamplingRate	48
4.3.2.28	RC MaxSamplingRate	49
4.3.3	Read Setup Commands	49
4.3.3.1	RS Accessory	50
4.3.3.2	RS Filter	50
4.3.3.3	RS Aperture	50
4.3.3.4	RS Mode	51
4.3.3.5	RS RangeMode	51
4.3.3.6	RS Range	51
4.3.3.7	RS ExposureMode	51
4.3.3.8	RS Exposure	52
4.3.3.9	RS SyncMode	52
4.3.3.10	RS SyncFreq	52
4.3.3.11	RS ExposureX	53
4.3.3.12	RS MatrixMode	53
4.3.3.13	RS UserCalibMode	53
4.3.3.14	RS Matrix	53
4.3.3.15	RS Match	54
4.3.3.16	RS Speed	54
4.3.3.17	RS SamplingRate	54
4.3.3.18	RS MaxFreqFlickerSearch	55
4.3.3.19	RS CMF	55
4.3.4	Measure Commands	55
4.3.4.1	M	55
4.3.4.2	MT	56
4.3.4.3	MA	56
4.3.4.4	MF	56
4.4	Configure Commands	56
4.4.1	Configure Calibration Commands	57
4.4.1.1	CC Matrix	57
4.4.1.2	CC Match	57
4.5	Response Codes	58

<b>CONTENTS</b>	<b>1</b>
<b>5   Deprecated List</b>	<b>61</b>
<b>Index</b>	<b>61</b>

---



# Chapter 1

## Introduction

The following sections forms the Remote Communications user manual:

- Section [Remote Communication](#) discusses how to download install drivers for your platform.
- Section [Getting Started](#) tells you how to generate your first piece of communication quickly.
- Section [Remote Commands](#) lists all the available commands and it's uses.





## Chapter 2

# Remote Communication

The CRI colorimeters comes standard with a well documented, easy to learn, pseudo english language based, command interpreter to control all aspects of it's operation using a computer, tablet or a smart device, making it a easy for customers to create their own software dedicated to perform specific measurement tasks or for inclusion in an Automated Test Environment. In addition, a fully documented communication and calculation libraries with numerous real world sample templates are included as starting points for building your own software tools using any of the modern computer development environments such as Visual Studio, Xcode etc.

### 2.1 Installing Drivers

On some platforms a device driver is necessary for the instrument to be recognized by the operating system. This will be available as part of your installer.

#### 2.1.1 Windows

To identify the Serial Port assigned,



1. Click on the Start Orb or button.

2. In the Start Search box type: device manager and then press enter and you should see something similar to the below example.




Figure 2.1: Device Manager

### 2.1.2 OSX

There is no USB drivers required for use with OSX

## 2.2 Remote Communication

A remote message consists of a *command* and an optional *key/value* followed by a [CF] or [LF] or [CF][LF] which is represented by a ← in this documentation. A command must begin with a root character. A root character also represents a command category.

#### Note

Remote messages are case-sensitive.

A format of a complete remote message is

```
[ROOT] [EXTENSION] [KEY] [VALUE] ←
```

There are four different command formats

- **Primary command**

A primary command consists of only the root character

Examples of primary commands are

---

E←

and

M←.

- **Secondary command**

Secondary command consists of root and extension characters.

Examples of a secondary command is

MF←

- **Tertiary command**

Tertiary commands consists of a secondary command followed by a **space** and a **key**.

Examples of tertiary commands are

>RC Time←

- **Complex command**

A complex command is a secondary command followed by a **key** and a **value** pair separated by **spaces**.

Example of complex command is

>SM Accessory 0←

The response of a remote message is detailed in [Response](#)

Detailed descriptions of remote messages explaining the *root*, *extension*, *key*, *value(s)* are described in the [Remote Commands](#) section.

Sample Command Syntax

>R\_\_ K, V1, V2←

The following table summarizes the high level categories in the Remote communication language.

Root	Description
E	<a href="#">Echo Command</a>
S	<a href="#">Setup Commands</a>
R	<a href="#">Read Commands</a>
M	<a href="#">Measure Commands</a>

### 2.2.1 Response

Every message is reciprocated with a success or a failure response.

#### Success

A successful message returns with the following response along with the message result if any.

OK: **code**: **response**←

Examples:

For the message

>SM Accessory 0←

your response would be as follows

OK:0:Accessory:No errors

## Failure

Error responses are in the following format

ER:*code:description:message*←

### Examples:

For the message

>SM Accessory←

if you typed incorrectly

>SM Accessory1← your response would be as follows

ER:-500:Invalid command:Accessory1←

## 2.2.2 Coding

The preferred method of dispatching commands to the instrument is sequential and polling to verify the command returns with a valid result or an error. Each command should return fairly quickly except for the measurement commands. For the measurement command the time out should be a multiple of the exposure time and the exposure multiplier/average. If an event driven approach is taken certain commands causes the command queue to ignore commands till the commands are processed. Care should be taken as a limited buffer is available to queue incoming commands. For eg. the M command ignores all other [M](#) command till it finishes the current measurement. Also another example is the [RM Spectrum](#) command, here a small delay is required preferable 200ms to allow the command to finish processing before a new command is queued.

---

## Chapter 3

# Getting Started

### 3.1 Getting started with the Terminal Software

#### 3.1.1 Windows®

##### Windows PowerShell

```
C:\>powershell
PS> Get-WMIObject Win32_SerialPort | Select-Object Name,DeviceID,Description
```

##### Windows Management Instrumentation Command-line (WMIC)

```
C:\>wmic path Win32_SerialPort Where "Caption LIKE '%COM%'" Get DeviceID
```

##### 3.1.1.1 Hyper Terminal

Hyper Terminal

##### 3.1.1.2 Windows®: PuTTY

PuTTY

#### 3.1.2 OS X

##### 3.1.2.1 screen

screen

##### Note

While using screen you will have to close the screen session by using Ctrl + A + k or Ctrl + A + \ to disconnect free up any serial port resources.

OS X terminal, list the available serial devices

```
ls /dev/cu*
```

### 3.2 Getting started with Remote messaging

As soon as a connection is established in any of the terminal applications or using a programming language, the instrument is ready to accept commands.

If you are running in an interactive shell it's best if you turn on duplex mode using the [Echo Command](#). This command will prompt >.

The remote command

>M

Takes a measurement.

See Also

[M](#)

>RM xy

Reports a measurement's CIE 1931 xy parameters.

See Also

[Read Measurement Commands](#)

---

## Chapter 4

# Remote Commands

### 4.1 Echo Command

Since

Firmware 1.03.

When typing out commands on a terminal sometimes it's hard to know what you type. To make communication intuitive the instrument can be instructed to remotely echo the sent commands received back to the sender. The E command toggles the ECHO of characters from the remote instrument.

To turn on/off ECHO send the E character

E←

Result: ECHO is **ON** the > prompt is echoed on to the terminal. ECHO is **OFF** nothing is echoed on the terminal.

Warning

When communicating with an automated system ensure that the ECHO is turned OFF in order to prevent confusing the controlling software.

### 4.2 Setup Commands

The setup command set is used to specify measurement and hardware (instrument) properties.

See Also

[Read Commands.](#)

S Commands	Description
SM	<a href="#">Measurement Setup Commands</a>
SH	<a href="#">Hardware Setup Commands</a>
SC	<a href="#">Computation Setup Commands</a>

#### 4.2.1 Measurement Setup Commands

The measurement setup command set is used to specify measurement properties.

**Quick Reference**



Key	Description
Accessory	Selects the accessory using its identifier
Filter1	Selects the 1 <sup>st</sup> filter using its identifier
Filter2	Selects the 2 <sup>nd</sup> filter using its identifier
Filter3	Selects the 3 <sup>rd</sup> filter using its identifier
Aperture	Selects the aperture using its identifier
Mode	Selects the instrument mode
ExposureMode	Selects the exposure mode
Exposure	Selects the exposure
MaxAutoExposure	Selects the maximum auto exposure limit
RangeMode	Selects the range mode
Range	Selects the range
SyncMode	Selects the sync mode
SyncFreq	Selects the Sync Frequency
ExposureX	Selects the exposure multiplier
MatrixMode	Selects the matrix mode. <i>Deprecated as of Firmware 1.16.</i>
UserCalibMode	Selects the user calibration mode.
Matrix	Selects the calibration matrix
Match	Selects the match calibration
Speed	Selects the speed
SamplingRate	Selects the sampling rate
MaxFreqFlickerSearch	Sets the Maximum Frequency Flicker Search
CMF	Selects the color matching function
Reset	Resets the instrument setup parameters to factory defaults.

#### 4.2.1.1 SM Accessory

##### Syntax

```
SM Accessory [ID]←
```

Selects the accessory using its ID.

##### See Also

[RC Accessory](#) lists all available accessories and their IDs.

##### Since

Firmware 1.04.

##### Examples

```
>SM Accessory 0
OK:0:Accessory:No errors
>SM Accessory -1
ER:-506:Accessory:Index doesn't select an Accessory
```

#### 4.2.1.2 SM Filter1

##### Syntax

```
SM Filter1 [ID]←
```

Selects the 1<sup>st</sup> filter using its ID.

**See Also**

[RC Filter](#) lists all available filters and their IDs.

**Since**

Firmware 1.04.

**Examples**

```
>SM Filter1 3
OK:0:Filter1:No errors
>SM Filter1 0
ER:-507:Filter1:Index doesn't select a Filter
```

**4.2.1.3 SM Filter2****Syntax**

```
SM Filter2 [ID]←
```

Selects the 2<sup>nd</sup> filter using its ID.

**See Also**

[RC Filter](#) lists all available filters and their IDs.

**Since**

Firmware 1.04.

**Examples**

```
>SM Filter2 3
OK:0:Filter2:No errors
```

**4.2.1.4 SM Filter3****Syntax**

```
SM Filter3 [ID]←
```

Selects the 3<sup>rd</sup> filter using its ID.

**See Also**

[RC Filter](#) lists all available filters and their IDs.

**Since**

Firmware 1.04.

**Examples**

```
>SM Filter3 3
OK:0:Filter3:No errors
```

---

#### 4.2.1.5 SM Aperture

##### Syntax

```
SM Aperture [ID]←
```

Selects the aperture using its ID.

##### See Also

[RC Aperture](#) lists all available apertures and their IDs.

##### Since

Firmware 1.04.

##### Examples

```
>SM Aperture 0
OK:0:SM Aperture:No errors
>SM Aperture -1
ER:-554:SM Aperture:Invalid argument:-1
>SM Aperture 1
ER:-515:SM Aperture:Index doesn't select an Aperture
```

#### 4.2.1.6 SM Mode

##### Syntax

```
SM Mode [ID]←
```

Selects the Instrument Mode using its ID.

##### See Also

[RC Mode](#) lists all available instrument modes and their IDs.

##### Since

Firmware 1.16.

##### Examples

```
>SM Mode 0
OK:0:SM Mode:No errors
>SM Mode -1
ER:-560:SM Mode:Invalid Instrument Mode
```

#### 4.2.1.7 SM ExposureMode

##### Syntax

```
SM ExposureMode [ID]←
```

Selects the exposure mode using its ID.

##### See Also

[RC ExposureMode](#) lists all available exposure modes and their IDs.

---

**Since**

Firmware 1.04.

**Examples**

```
>SM ExposureMode 0
OK:0:ExposureMode:No errors
>SM ExposureMode -1
ER:-518:ExposureMode:Invalid Exposure Mode
```

**4.2.1.8 SM Exposure****Syntax**

```
SM Exposure [exposure in msecs]←
```

Selects the exposure in msecs. This exposure is used only when the Exposure Mode is set to manual.

**See Also**

[RC MinExposure](#) and [RC MaxExposure](#) provides the range of exposures.

**Since**

Firmware 1.04.

**Examples**

```
>SM Exposure 10
OK:0:Exposure:No errors
>SM Exposure 1000000
ER:-519:Exposure:Invalid Exposure value
```

**4.2.1.9 SM MaxAutoExposure****Syntax**

```
SM MaxAutoExposure [exposure in msecs]←
```

Selects the maximum auto exposure limit in msecs. The exposure is limited to the maximum auto exposure limit when the Exposure Mode is set to auto.

**See Also**

[RC MinExposure](#) and [RC MaxExposure](#) provides the range of max auto exposure limit.

**Since**

Firmware 1.26.

**Examples**

```
>SM MaxAutoExposure 400
OK:0:MaxAutoExposure:No errors
```

---

#### 4.2.1.10 SM RangeMode

##### Syntax

```
SM RangeMode [ID]←
```

Selects the range mode using its ID.

##### See Also

[RC RangeMode](#) lists all available range modes and their IDs.

##### Since

Firmware 1.04.

##### Examples

```
>SM RangeMode 0
OK:0:RangeMode:No errors
>SM RangeMode -1
ER:-512:RangeMode:Invalid Range mode
```

#### 4.2.1.11 SM Range

##### Syntax

```
SM Range [ID]←
```

Selects the range using its ID. This is used only when the Range Mode is set to manual.

##### See Also

[RC Range](#) lists all available ranges and their IDs.

##### Since

Firmware 1.04.

##### Examples

```
>SM Range 1
OK:0:Range:No errors
>SM Range -1
ER:-513:Range:Invalid Range index
```

#### 4.2.1.12 SM SyncMode

##### Syntax

```
SM SyncMode [ID]←
```

Selects the sync mode using its ID.

##### See Also

[RC SyncMode](#) lists all available sync modes and their IDs.

---

**Since**

Firmware 1.04.

**Examples**

```
>SM SyncMode 0
OK:0:SyncMode:No errors
>SM SyncMode -1
ER:-521:SyncMode:Invalid Sync Mode
```

**4.2.1.13 SM SyncFreq****Syntax**

```
SM SyncMode [frequency in Hz.]←
```

Selects the Sync Frequency in Hz. This is used only when the Sync Mode is set to manual.

**Since**

Firmware 1.04.

**Examples**

```
>SM SyncFreq 10
OK:0:SyncFreq:No errors
>SM SyncFreq 0
ER:-522:SyncFreq:Invalid User Sync Frequency
```

**4.2.1.14 SM ExposureX****Syntax**

```
SM ExposureX [exposure multiplier]←
```

Selects the exposure multiplier.

**See Also**

[RC MinExposureX](#) and [RC MaxExposureX](#) provides the range of values that can be set.

**Since**

Firmware 1.04.

**Examples**

```
>SM ExposureX 1
OK:0:ExposureX:No errors
>SM ExposureX 0
ER:-514:ExposureX:Invalid Exposure Multiplier
```

**4.2.1.15 SM MatrixMode****Syntax**

```
SM MatrixMode [ID]←
```

Selects the matrix mode using its ID.

---

**Matrix Mode values**

- 0: Disabled
- 1: Enabled

**Since**

Firmware 1.04.

**Deprecated** [SM MatrixMode](#) has been deprecated for [SM UserCalibMode](#) as of Firmware 1.16

**Examples**

```
>SM MatrixMode 0
OK:0:MatrixMode:No errors
>SM MatrixMode -1
ER:-552:MatrixMode:Invalid Matrix Mode
```

**4.2.1.16 SM UserCalibMode****Syntax**

```
SM UserCalibMode [ID]←
```

Selects the user calibration mode using its ID.

**User Calibration Mode values**

- 0: None or Disabled
- 1: Matrix Calibration
- 2: Match Calibration

**Since**

Firmware 1.16.

**Examples**

```
>SM UserCalibMode 0
OK:0:UserCalibMode:No errors
>SM UserCalibMode -1
ER:-552:SM UserCalibMode:Invalid User Calibration Mode
```

**4.2.1.17 SM Matrix****Syntax**

```
SM Matrix [ID]←
```

Selects the calibration matrix by its ID.

**See Also**

[RC MatrixCalib](#) lists all available matrices and their IDs.

---

**Since**

Firmware 1.04.

**Examples**

```
>SM Matrix 0
OK:0:Matrix:No errors
>SM Matrix -1
ER:-553:SM Matrix:Invalid Matrix ID
```

**4.2.1.18 SM Match****Syntax**

```
SM Match [ID]←
```

Selects the Match Calibration set to be used by its ID.

**See Also**

[RC MatchCalib](#) lists all available match calibrations and their IDs.

**Since**

Firmware 1.16.

**Examples**

```
>SM Match 0
OK:0:Match:No errors
>SM Match -1
ER:-557:SM Match:Invalid Match ID
```

**4.2.1.19 SM Speed****Syntax**

```
SM Speed [ID]←
```

Selects the Speed to be used by its ID.

**See Also**

[RC Speed](#) lists all available speeds and their IDs.

**Note**

This command is only valid if the [RC InstrumentType](#) is 2 (Spectroradiometer).

**Since**

Firmware 1.17.

**Examples**

```
>SM Speed 0
OK:0:Speed:No errors
>SM Speed -1
ER:-557:SM Speed:Invalid Speed ID
```

---



#### 4.2.1.20 SM SamplingRate

##### Syntax

```
SM SamplingRate [frequency in Hz.]←
```

Selects the Sampling Rate in Hz. This is used only when the Instrument Mode is set to Flicker or Response Time.

##### Since

Firmware 1.19

##### Examples

```
>SM SamplingRate 220
OK:0:SamplingRate:No errors
>SM SamplingRate 0
ER:-522:SamplingRate:Invalid Sampling Rate
```

#### 4.2.1.21 SM MaxFreqFlickerSearch

##### Syntax

```
SM MaxFreqFlickerSearch [frequency in Hz.]←
```

Selects the Maximum Frequency Flicker Search in Hz.

##### Since

Firmware 1.19

##### Examples

```
>SM MaxFreqFlickerSearch 220
OK:0:SM MaxFreqFlickerSearch:No errors
>SM MaxFreqFlickerSearch -1
ER:-524:SM MaxFreqFlickerSearch:Invalid MaxFreqFlickerSearch
```

#### 4.2.1.22 SM CMF

##### Syntax

```
SM CMF [index]←
```

Selects the Color Matching Function to be used for the measurement calculations.

##### Note

This command is only valid if the [RC InstrumentType](#) is 2 (Spectroradiometer).

##### Since

Firmware 1.26

##### Examples

```
>SM CMF 1
OK:0:CMF:No errors
```

---

## 4.2.1.23 SM Reset

v

## Syntax

SM Reset← Resets the instrument setup parameters to factory defaults.

## Since

Firmware 1.36

## Examples

```
>SM Reset
OK:0:SM Reset:No errors
```

For the CR-100, the measurement Setup Reset command "SM Reset" is equivalent to the following SM commands-:

Command	Reset value
SM Mode	0 = Colorimeter
SM Accessory	0 = Installed objective lens
SM Filter1	0 = No filter selected
SM Filter2	0 = No filter selected
SM Filter3	0 = No filter selected
SM SyncMode	0 = No synchronization
SM SyncFreq	60 Hz
SM RangeMode	0 = Automatic
SM Range	0 = Range A
SM Exposure	1 = 1 millisecond exposure
SM MaxAutoExposure	500 = 500 milliseconds
SM ExposureX	1
SM UserCalibMode	0 = None
SM SamplingRate	1000 Hz
SM MaxFreqFlickerSearch	120 Hz

## 4.2.2 Hardware Setup Commands

The hardware setup command set is used to specify instrument properties. This section is reserved for future expansion.

## Quick Reference

Key	Description
-----	-------------

## 4.2.3 Computation Setup Commands

The computation setup command set is used to specify computation properties. This section is reserved for future expansion.

## Quick Reference

Key	Description
CMF	Selects the color matching function

## 4.2.3.1 SC CMF

## Syntax

SC CMF [index]←

Selects the Color Matching Function table index.

#### Note

This command is only valid if the [RC InstrumentType](#) is 2 (Spectroradiometer).  
This command also recalculates and modifies the current measurement.

CMF indices are as follows

- **0**: default CIE 1931.
- **1**: user CMF 1.
- **2**: user CMF 2.
- **3**: user CMF 3.

#### Since

Firmware 1.26

#### Examples

```
>SC CMF 0
OK:0:SC CMF:No errors
```

## 4.3 Read Commands

The read command set is used to retrieve properties from the instrument.

R Commands	Description
RM	<a href="#">Read Measurement Commands</a>
RS	<a href="#">Read Setup Commands</a>
RC	<a href="#">Read Configuration Commands</a>

### 4.3.1 Read Measurement Commands

The read measurement command set is used to retrieve properties from the last measurement.

#### Warning

The following commands are only valid after a measurement using a [M](#) or [MT](#)

#### Quick Reference

Key	Description
ID	Instrument Identifier (Serial Number)
Model	Instrument Model
Time	Date and Time stamp of measurement
Accessory	Accessories used during capture
Filter	Filters used during capture
Aperture	Aperture used during capture
Mode	Instrument Mode
ExposureMode	Exposure Mode

Exposure	Exposure time in milliseconds
MaxAutoExposure	Maximum auto exposure limit in milliseconds
RangeMode	Range Mode
Range	Range used during capture
SyncMode	Sync Mode during capture
SyncFreq	Sync Frequency
ExposureX	Exposure Multiplier
MatrixMode	Calibration matrix enabled/disabled. <i>Deprecated as of Firmware 1.16.</i>
UserCalibMode	User Calibration mode selection.
Matrix	Current calibration matrix ID
Match	Current match calibration ID
Speed	Speed used during capture
X	2°Tristimulus X
X10	10°Tristimulus X
Y	2°Tristimulus Y
Y10	10°Tristimulus Y
Z	2°Tristimulus Z
Z10	10°Tristimulus Z
XYZ	2°Tristimulus data
XYZ10	10°Tristimulus data
xy	CIE 1931 xy data (2°)
xy10	CIE 1964 xy data (10°)
uv	CIE 1960 uv data
upvp	CIE 1976 u' v' data
CCT	Correlated color temperature (CCT) in °Kelvin and uv Deviation from the planckian locus.
Warnings	Measurement warnings
Ye	
Yv	
Radiometric	
Spectrum	Spectral response
Speed	Speed settings
Temporal	Temporal response
TemporalY	Corrected temporal response
SamplingRate	Sampling rate
CMF	Color matching function used

#### 4.3.1.1 RM ID

##### Syntax

```
RM ID←
```

Retrieves the instrument identifier (Serial Number) stored during the last capture.

##### Since

Firmware 1.04.

##### Examples

```
>RM ID
OK:0:RM ID:A00102
```

#### 4.3.1.2 RM Model

##### Syntax

```
RM Model←
```

Retrieves the instrument model stored during the last capture.

##### Since

Firmware 1.04.

##### Examples

```
>RM Model  
OK:0:RM Model:CR-100
```

#### 4.3.1.3 RM Time

##### Syntax

```
RM Time←
```

Retrieves the date and time stamp of measurement

##### Since

Firmware 1.04.

##### Examples

```
>RM Time  
OK:0:RM Time:NA
```

#### 4.3.1.4 RM Accessory

##### Syntax

```
RM Accessory←
```

Retrieves the name of the accessory used during the last capture.

##### See Also

[RC Filter](#) to get the accessory attributes.

##### Since

Firmware 1.04.

##### Examples

```
>RM Accessory  
OK:0:RM Accessory:Standard
```

#### 4.3.1.5 RM Filter

##### Syntax

```
RM Filter←
```

Retrieves the names of the filters used during the last capture, as a comma separated values.

---

**See Also**

[RC Filter](#) to get the filter attributes.

**Since**

Firmware 1.04.

**Examples**

```
>RM Filter
OK:0:RM Filter:None
>RM Filter
OK:0:RM Filter:ND-100-1,ND-100-3
```

**4.3.1.6 RM Aperture****Syntax**

```
RM Aperture←←
```

Retrieves the name of the aperture used during last capture.

**See Also**

[RC Aperture](#) to get the aperture attributes.

**Since**

Firmware 1.04.

**Examples**

```
>RM Aperture
OK:0:RM Aperture:5 deg
```

**4.3.1.7 RM Mode****Syntax**

```
RM Mode←←
```

Retrieves the name of the instrument mode used during last capture.

**Since**

Firmware 1.16.

**Examples**

```
>RM Mode
OK:0:RM Mode:Colorimeter
```

**4.3.1.8 RM ExposureMode****Syntax**

```
RM ExposureMode←←
```

Retrieves the name of the exposure mode used during last capture.

---

**See Also**

[RC ExposureMode](#) to get the exposure mode attributes.

**Since**

Firmware 1.04.

**Examples**

```
>RM ExposureMode
OK:0:RM ExposureMode:Auto
```

**4.3.1.9 RM Exposure****Syntax**

```
RM Exposure←
```

Exposure time in milliseconds

**Since**

Firmware 1.04.

**Examples**

```
>RM Exposure
OK:0:RM Exposure:111.622 msec
```

**4.3.1.10 RM MaxAutoExposure****Syntax**

```
RM MaxAutoExposure←
```

Maximum auto exposure limit in milliseconds

**Since**

Firmware 1.26.

**Examples**

```
>RM MaxAutoExposure
OK:0:RM MaxAutoExposure:449.999 msec
```

**4.3.1.11 RM RangeMode****Syntax**

```
RM RangeMode←
```

Retrieves the name of the range mode used during last capture.

**See Also**

[RC RangeMode](#) to get the range mode attributes.

---

**Since**

Firmware 1.04.

**Examples**

```
>RM RangeMode
OK:0:RM RangeMode:Auto
```

**4.3.1.12 RM Range****Syntax**

```
RM Range←
```

Retrieves the name of the range used during last capture.

**See Also**

[RC Range](#) to get the range attributes.

**Since**

Firmware 1.04.

**Examples**

```
>RM Range
OK:0:RM Range:D
```

**4.3.1.13 RM SyncMode****Syntax**

```
RM SyncMode←
```

Retrieves the name of the sync mode used during last capture.

**See Also**

[RC SyncMode](#) to get the sync mode attributes.

**Since**

Firmware 1.04.

**Examples**

```
>RM SyncMode
OK:0:RM SyncMode:None
```

**4.3.1.14 RM SyncFreq****Syntax**

```
RM SyncFreq←
```

Retrieves the sync frequency during last capture.

---



**Warning**

This command should not be used if the [RM SyncMode](#) command returns **None**.

**Since**

Firmware 1.04.

**Examples**

```
>RM SyncFreq
OK:0:RM SyncFreq:0.00 Hz
```

**4.3.1.15 RM ExposureX****Syntax**

```
RM ExposureX←
```

Retrieves the exposure multiplier used during last capture.

**Since**

Firmware 1.04.

**Examples**

```
>RM ExposureX
OK:0:RM ExposureX:1
```

**4.3.1.16 RM MatrixMode****Syntax**

```
RM MatrixMode←
```

Retrieves the matrix mode used during last capture. If matrix mode is enabled then a calibration matrix was applied the data. [RM Matrix](#) will retrieve the calibration matrix ID.

**Since**

Firmware 1.04.

**Deprecated** [RM MatrixMode](#) has been deprecated for [RM UserCalibMode](#) as of Firmware 1.16

**Examples**

```
>RM MatrixMode
OK:0:RM MatrixMode:Disabled
```

**4.3.1.17 RM UserCalibMode****Syntax**

```
RM UserCalibMode←
```

Reports the User Calibration Mode that was in effect when the last measurement was taken.

If user calibration mode is set to **Matrix** then a calibration matrix was applied the data. [RM Matrix](#) will retrieve the calibration matrix ID. If user calibration mode is set to **Match** then a match calibration was applied the data. [RM Match](#) will retrieve the match calibration ID.

---

**Since**

Firmware 1.16.

**Examples**

```
>RM UserCalibMode
OK:0:RM UserCalibMode:None
```

**4.3.1.18 RM Matrix****Syntax**

```
RM Matrix←
```

Retrieves the matrix calibration ID used during last capture.

**Warning**

This command should be used if the [RM UserCalibMode](#) is **Match**.

**Since**

Firmware 1.04.

**Examples**

```
>RM Matrix
OK:0:RM Matrix:N
>RM Matrix
OK:0:RM Matrix:1
```

**4.3.1.19 RM Match****Syntax**

```
RM Match←
```

Reports the Match Calibration set ID that was in effect when the last measurement was taken.

**Warning**

This command should be used if the [RM UserCalibMode](#) is **Match**.

**Since**

Firmware 1.16.

**Examples**

```
>RM Match
OK:0:RM Match:0
>RM Match
OK:0:RM Match:1
```

**4.3.1.20 RM Speed****Syntax**

```
RM Speed←
```

Retrieves the name of the speed used during last capture.

---

**See Also**

[RC Speed](#) to get the speed attributes.

**Note**

This command is only valid if the [RC InstrumentType](#) is 2 (Spectroradiometer).

**Since**

Firmware 1.17.

**Examples**

```
>RM Speed
OK:0:RM Speed:Normal
```

**4.3.1.21 RM X****Syntax**

```
RM X←
```

Retrieves the (2°) tristimulus X data of the last capture.

**Since**

Firmware 1.04.

**Examples**

```
>RM X
OK:0:RM X:1.737e+00
```

**4.3.1.22 RM X10****Syntax**

```
RM X10←
```

Retrieves the (10°) tristimulus X data of the last capture.

**Note**

This command is only valid if the [RC InstrumentType](#) is 2 (Spectroradiometer).

**Since**

Firmware 1.18.

**Examples**

```
>RM X10
OK:0:RM X10:1.737e+00
```

---

#### 4.3.1.23 RM Y

##### Syntax

RM Y←

Retrieves the (2°) tristimulus Y data of the last capture.

##### Since

Firmware 1.04.

##### Examples

```
>RM Y
OK:0:RM Y:1.685e+00
```

#### 4.3.1.24 RM Y10

##### Syntax

RM Y10←

Retrieves the (10°) tristimulus Y data of the last capture.

##### Note

This command is only valid if the [RC InstrumentType](#) is 2 (Spectroradiometer).

##### Since

Firmware 1.18.

##### Examples

```
>RM Y10
OK:0:RM Y10:1.685e+00
```

#### 4.3.1.25 RM Z

##### Syntax

RM Z←

Retrieves the (2°) tristimulus Z data of the last capture.

##### Since

Firmware 1.04.

##### Examples

```
>RM Z
OK:0:RM Z:1.830e+00
```

#### 4.3.1.26 RM Z10

##### Syntax

RM Z10←

Retrieves the (10°) tristimulus Z data of the last capture.

---

**Note**

This command is only valid if the [RC InstrumentType](#) is 2 (Spectroradiometer).

**Since**

Firmware 1.18.

**Examples**

```
>RM Z10
OK:0:RM Z10:1.830e+00
```

**4.3.1.27 RM XYZ****Syntax**

```
RM XYZ←
```

Retrieves the (2°) tristimulus XYZ data of the last capture as comma separated values.

**Since**

Firmware 1.04.

**Examples**

```
>RM XYZ
OK:0:RM XYZ:1.737e+00,1.685e+00,1.830e+00
```

**4.3.1.28 RM XYZ10****Syntax**

```
RM XYZ10←
```

Retrieves the (10°) tristimulus XYZ data of the last capture as comma separated values.

**Note**

This command is only valid if the [RC InstrumentType](#) is 2 (Spectroradiometer).

**Since**

Firmware 1.18.

**Examples**

```
>RM XYZ10
OK:0:RM XYZ10:1.737e+00,1.685e+00,1.830e+00
```

**4.3.1.29 RM xy****Syntax**

```
RM xy←
```

Retrieves the CIE 1931 xy data (2°) of the last capture as comma separated values.

---

**Since**

Firmware 1.04.

**Examples**

```
>RM xy
OK:0:RM xy:0.3308,0.3208
```

**4.3.1.30 RM xy10****Syntax**

```
RM xy10←
```

Retrieves the CIE 1964 xy data (10°) of the last capture as comma separated values.

**Note**

This command is only valid if the [RC InstrumentType](#) is 2 (Spectroradiometer).

**Since**

Firmware 1.18.

**Examples**

```
>RM xy
OK:0:RM xy:0.3308,0.3208
```

**4.3.1.31 RM uv****Syntax**

```
RM uv←
```

Retrieves the CIE 1960 uv data of the last capture as comma separated values.

**Since**

Firmware 1.04.

**Examples**

```
>RM uv
OK:0:RM uv:0.2138,0.3110
```

**4.3.1.32 RM upvp****Syntax**

```
RM upvp←
```

Retrieves the CIE 1976 u' v' data of the last capture as comma separated values.

**Since**

Firmware 1.04.

**Examples**

```
>RM upvp
OK:0:RM upvp:0.2138,0.4666
```

---

#### 4.3.1.33 RM CCT

##### Syntax

```
RM CCT←
```

Retrieves the correlated color temperature (CCT) in °Kelvin and the uv deviation from the planck's locus of the last capture as comma separated values.

##### Since

Firmware 1.04.

##### Examples

```
>RM CCT  
OK:0:RM CCT:5577,-0.0100
```

#### 4.3.1.34 RM Warnings

##### Syntax

```
RM Warnings←
```

##### Since

Firmware 1.04.

Retrieves the measurement warnings reported during the last capture.

##### Examples

```
>RM Warnings  
OK:0:RM Warnings:0
```

#### 4.3.1.35 RM Yv

##### Syntax

```
RM Yv←
```

##### Since

Firmware 1.17.

Retrieves the scotopic Luminance.

##### Examples

```
>RM Yv  
OK:0:RM Yv:0
```

#### 4.3.1.36 RM Radiometric

##### Syntax

```
RM Radiometric←
```

---

**Since**

Firmware 1.17.

Retrieves the Units of radiometric data, Radiometric power and Photon Radiometric power.

values for radiometric types are

- **0**: Radiance.
- **1**: Irradiance.
- **2**: Radiant Intensity.
- **3**: Radiant Flux.

**Examples**

```
>RM Radiometric
OK:0:RM Radiometric:0,3.209e-01,8.835e+17
```

**4.3.1.37 RM Spectrum****Syntax**

```
RM Spectrum←
```

**Since**

Firmware 1.17.

Retrieves the spectral measurement. If a valid spectral reading is present the first line will indicate start, end, delta wavelengths and the number of points. Followed by a spectral point on each line.

**Note**

After issuing this command a small delay (200ms) is required before another command is issued. The preferred method of commands is polling to see if the command returns but in case an event driven approach is taken this delay is required.

**Examples**

```
>RM Spectrum
OK:0:RM Spectrum:380.0,780.0,2.0,201
2.119e-24
1.913e-24
.
.
.
```

**4.3.1.38 RM Temporal****Syntax**

```
RM Temporal←
```

**Since**

Firmware 1.19.

Retrieves the temporal measurement. If a valid temporal reading is present the first line will indicate sampling rate and the number of points. Followed by a temporal point on each line.

---



**Note**

After issuing this command a small delay (200ms) is required before another command is issued. The preferred method of commands is polling to see if the command returns but in case an event driven approach is taken this delay is required.

**Examples**

```
>RM Temporal
OK:0:RM Temporal:220.0,1024
20827
20841
20851
.
.
.
```

**4.3.1.39 RM TemporalY****Syntax**

```
RM TemporalY←
```

**Since**

Firmware 1.20.

Retrieves the calibrated temporal measurement. If a valid temporal reading is present the first line will indicate sampling rate and the number of points. Followed by a calibrated temporal point on each line.

**Note**

After issuing this command a small delay (200ms) is required before another command is issued. The preferred method of commands is polling to see if the command returns but in case an event driven approach is taken this delay is required.

**Examples**

```
>RM TemporalY
OK:0:RM TemporalY:220.0,1024
1.838e+00
1.826e+00
1.855e+00
.
.
.
```

**4.3.1.40 RM SamplingRate****Syntax**

```
RM SamplingRate←
```

**Since**

Firmware 1.19.

Retrieves the sampling rate used during the temporal measurement.

---

**Note**

After issuing this command a small delay (200ms) is required before another command is issued. The preferred method of commands is polling to see if the command returns but in case an event driven approach is taken this delay is required.

**Examples**

```
>RM SamplingRate
OK:0:RM SamplingRate:200.0
```

**4.3.1.41 RM CMF****Syntax**

```
RM CMFv←
```

**Since**

Firmware 1.26.

Retrieves the Color Matching Function used for calculations.

**Note**

This command is only valid if the [RC InstrumentType](#) is 2 (Spectroradiometer).

**Examples**

```
>RM CMF
OK:0:RM CMF:0
```

**4.3.2 Read Configuration Commands****Quick Reference**

Key	Description
ID	Instrument Identifier/Serial Number
Model	Instrument Model Number
InstrumentType	Instrument type
Accessory	Available accessories
Filter	Available filters
Aperture	Available apertures
ExposureMode	Available exposure modes
RangeMode	Available range modes
Range	Available ranges
SyncMode	Available sync modes
Firmware	Firmware version
MatrixMode	Available matrix modes. <i>Deprecated as of Firmware 1.16.</i>
Matrix	Available matrix calibrations
MatrixCalibration	Available matrix calibration factors
MinExposure	Minimum exposure time

MaxExposure	Maximum exposure time
MinSyncFreq	Minimum user selectable sync frequency
MaxSyncFreq	Maximum user selectable sync frequency
MinExposureX	Minimum exposure multiplier
MaxExposureX	Maximum exposure multiplier
Speed	Speed
MinSamplingRate	Minimum sampling rate
MaxSamplingRate	Maximum sampling rate

#### 4.3.2.1 RC ID

##### Syntax

```
RC ID←
```

Reads the instrument Identifier (Serial Number) from the configuration.

##### Since

Firmware 1.04.

##### Examples

```
>RC ID
OK:0:RC ID:A00102
```

#### 4.3.2.2 RC Model

##### Syntax

```
RC Model←
```

Reads the instrument Model from the configuration.

##### Since

Firmware 1.04.

##### Examples

```
>RC Model
OK:0:RC Model:CR-100
```

#### 4.3.2.3 RC InstrumentType

##### Syntax

```
RC InstrumentType←
```

Reads the instrument type from the configuration.

values for instrument types are

- **0**: Photometer.
- **1**: Colorimeter.
- **2**: Spectroradiometer.

**Since**

Firmware 1.17.

**Examples**

```
>RC InstrumentType
OK:0:RC InstrumentType:2
```

**4.3.2.4 RC Accessory****Syntax**

```
RC Accessory←
```

List the available accessories that are configured.

The response returns the number of accessories in the last parameter (4<sup>th</sup>) of the first line. Every line after that lists the accessory ID, name and type.

**Accessory types**

- **Radiance:** Radiance
- **Irradiance:** Irradiance
- **Rad. Flux:** Radiant Flux
- **Rad. Intensity:** Radiant Intensity
- **NA:** Not applicable
- **Unassigned:** Unassigned

**Since**

Firmware 1.04.

**Examples**

```
>RC Accessory
OK:0:RC Accessory:3
0,Standard,Radiance
1,IR-100,Irradiance
2,IS-101,Rad. Flux
```

**4.3.2.5 RC Filter****Syntax**

```
RC Filter←
```

Lists the available filters that are configured.

The response returns the number of filters in the last parameter(4<sup>th</sup>) of the first line. Every line after that lists the filter ID, name and type.

**Filter types**

- **Radiance:** Radiance
  - **Irradiance:** Irradiance
  - **Rad. Flux:** Radiant Flux
  - **Rad. Intensity:** Radiant Intensity
  - **NA:** Not applicable
  - **Unassigned:** Unassigned
-

**Since**

Firmware 1.04.

**Examples**

```
>RC Filter
OK:0:RC Filter:5
3,ND-100-1,Radiance
4,ND-100-2,Radiance
5,ND-100-3,Radiance
6,ND-100-0.3,Radiance
7,ND-100-0.7,Radiance
```

**4.3.2.6 RC Aperture****Syntax**

RC Aperture←

Lists the available apertures that are configured.

The response returns the number of apertures in the last parameter(4<sup>th</sup>) of the first line. Every line after that lists the aperture ID and name.

**Since**

Firmware 1.04.

**Examples**

```
>RC Aperture
OK:0:RC Aperture:1
0,5 deg
```

**4.3.2.7 RC Mode****Syntax**

RC Mode←

List the available Instrument Modes available.

The response returns the number of instruments modes in the last parameter(4<sup>th</sup>) of the first line. Every line after that lists the instrument mode ID and name.

**Since**

Firmware 1.16.

**Exposure Modes**

- Colorimeter
- Flicker
- Response Time

**Examples**

```
>RC Mode
OK:0:RC Mode:3
0,Colorimeter
1,Flicker
2,Response Time
```

---

#### 4.3.2.8 RC ExposureMode

##### Syntax

```
RC ExposureMode←
```

List the available Exposure Modes configured.

The response returns the number of exposure modes in the last parameter(4<sup>th</sup>) of the first line. Every line after that lists the exposure mode ID and name.

##### Since

Firmware 1.04.

##### Exposure Modes

- Auto
- Fixed

##### Examples

```
>RC ExposureMode
OK:0:RC ExposureMode:2
0,Auto
1,Fixed
```

#### 4.3.2.9 RC RangeMode

##### Syntax

```
RC RangeMode←
```

Lists the available range modes configured.

The response returns the number of range modes in the last parameter(4<sup>th</sup>) of the first line. Every line after that lists the range mode ID and name.

##### Range Modes

- Auto
- Fixed

##### Since

Firmware 1.04.

##### Examples

```
>RC RangeMode
OK:0:RC RangeMode:2
0,Auto
1,Fixed
```

#### 4.3.2.10 RC Range

##### Syntax

```
RC Range←
```

---

Lists the available ranges configured.

The response returns the number of ranges in the last parameter(4<sup>th</sup>) of the first line. Every line after that lists the range ID and name.

#### Since

Firmware 1.04.

#### Examples

```
>RC Range
OK:0:RC Range:4
0,A
1,B
2,C
3,D
```

### 4.3.2.11 RC SyncMode

#### Syntax

```
RC SyncMode←
```

Lists the available Sync Modes configured.

The response returns the number of sync modes in the last parameter(4<sup>th</sup>) of the first line. Every line after that lists the sync mode ID and name.

#### SyncModes

- **None:** The instrument will not sync to the light source.
- **Auto:** The instrument syncs to the light source automatically.
- **Manual:** The instrument will sync to the frequency specified by [RS SyncFreq](#).

#### Since

Firmware 1.04.

Additional Sync modes to allow for easy selection based on the display technology

- **NTSC**
- **PAL**
- **CINEMA**

#### Since

Firmware 1.32

#### Examples

```
>RC SyncMode
OK:0:RC SyncMode:3
0,None
1,Auto
2,Manual
```

---

```
>RC SyncMode
OK:0:RC SyncMode:3
0,None
1,Auto
2,Manual
3,NTSC
4,PAL
5,CINEMA
```

#### 4.3.2.12 RC Firmware

##### Syntax

```
RC Firmware←
```

Reports the current Firmware version

##### Note

Use the firmware version to verify supported commands

##### Since

Firmware 1.04.

##### Examples

```
>RC Firmware
OK:0:RC Firmware:1.04
```

#### 4.3.2.13 RC MatrixMode

##### Syntax

```
RC MatrixMode←
```

Lists the available Matrix Modes configured.

The response returns the number of matrix modes in the last parameter(4<sup>th</sup>) of the first line. Every line after that lists the matrix mode ID and name.

##### MatrixModes

- **Disabled:** The instrument will not apply matrix calibration to the captured data.
- **Enabled:** The instrument will apply matrix calibration to the captured specified by [RS Matrix](#).

##### Since

Firmware 1.04.

**Deprecated** [RC MatrixMode](#) has been deprecated for [RC UserCalibMode](#) as of Firmware 1.16

##### Examples

```
>RC MatrixMode
OK:0:RC MatrixMode:2
0,Disabled
1,Enabled
```

---



#### 4.3.2.14 RC UserCalibMode

##### Syntax

```
RC UserCalibMode←
```

Lists the available user calibration modes configured for the current instrument.

The response returns the number of user calibration modes in the last parameter(4<sup>th</sup>) of the first line. Every line after that lists the user calibration mode ID and name.

##### User Calibration Modes

- **None:** The instrument will not apply matrix calibration to the captured data.
- **Matrix:** The instrument will apply matrix calibration to the captured specified by [RS Matrix](#).
- **Match:** The instrument will apply calibration match to the captured specified by [RS Match](#).

##### Since

Firmware 1.16.

##### Examples

```
>RC UserCalibMode
OK:0:RC UserCalibMode:3
0,None
1,Matrix
2,Match
```

#### 4.3.2.15 RC Matrix

##### Syntax

```
RC Matrix←
```

Lists the available matrix calibrations for the current selected accessory.

##### Since

Firmware 1.04.

##### Examples

```
>RC Matrix
OK:0:RC Matrix:None
>RC Matrix
OK:0:RC Matrix:1
0,Display Test
```

#### 4.3.2.16 RC Match

##### Syntax

```
RC Match←
```

Lists the available match calibration sets for the current instrument.

---

**Since**

Firmware 1.16.

**Examples**

```
>RC Match
OK:0:RC Match:None
>RC Match
OK:0:RC Match:1
0,Test
```

**4.3.2.17 RC MatrixCalibration****Syntax**

```
RC MatrixCalibration←
```

Lists the available Matrix Calibration sets.

**Since**

Firmware 1.04.

**Deprecated** [RC MatrixCalibration](#) has been deprecated for [RC MatrixCalib](#) as of Firmware 1.16

**Examples**

```
>RC MatrixCalibration
OK:0:RC MatrixCalibration:None
>RC MatrixCalibration
OK:0:RC MatrixCalibration:1
0,Display Test,1.030e+00,-1.363e-02,-8.051e-03,-2.175e-02,1.072e+00,1.-
203e-02,5.340e-02,3.940e-03,1.058e+00
```

**4.3.2.18 RC MatrixCalib****Syntax**

```
RC MatrixCalibration←
```

Lists the available Matrix Calibration sets.

**Since**

Firmware 1.16.

**Examples**

```
>RC MatrixCalib
OK:0:RC MatrixCalib:None
>RC MatrixCalib
OK:0:RC MatrixCalib:1
0,Display Test,1.030e+00,-1.363e-02,-8.051e-03,-2.175e-02,1.072e+00,1.-
203e-02,5.340e-02,3.940e-03,1.058e+00
```

**4.3.2.19 RC MatchCalib****Syntax**

```
RC MatchCalib←
```

Lists the available Match Calibration sets and factors.

---

**Since**

Firmware 1.16.

**Examples**

```
>RC MatchCalib
OK:0:RC MatchCalib:None
>RC MatchCalib
OK:0:RC MatchCalib:1
0,Test,5.292e-01,8.048e-01,7.837e-01
```

**4.3.2.20 RC MinExposure****Syntax**

```
RC MinExposure←
```

Lower limit of the exposure reported in milliseconds.

**See Also**

[RC MaxExposure](#) for upper limit for the exposure for Fixed exposure mode.

**Since**

Firmware 1.04.

**Examples**

```
>RC MinExposure
OK:0:RC MinExposure:1.0 msec
```

**4.3.2.21 RC MaxExposure****Syntax**

```
RC MaxExposure←
```

Upper limit of the exposure reported in milliseconds.

**See Also**

[RC MinExposure](#) for lower limit for the exposure for Fixed exposure mode.

**Since**

Firmware 1.04.

**Examples**

```
>RC MaxExposure
OK:0:RC MaxExposure:500.0 msec
```

**4.3.2.22 RC MinSyncFreq****Syntax**

```
RC MinSyncFreq←
```

Lower limit of the sync frequency reported in Hertz.

---

**See Also**

[RC MaxSyncFreq](#) for upper limit of the sync frequency for Manual sync mode.

**Since**

Firmware 1.04.

**Examples**

```
>RC MinSyncFreq
OK:0:RC MinSyncFreq:10.00 Hz
```

**4.3.2.23 RC MaxSyncFreq****Syntax**

```
RC MaxSyncFreq←
```

Upper limit of the sync frequency reported in Hertz.

**See Also**

[RC MinSyncFreq](#) for lower limit of the sync frequency for Manual sync mode.

**Since**

Firmware 1.04.

**Examples**

```
>RC MaxSyncFreq
OK:0:RC MaxSyncFreq:10000.00 Hz
```

**4.3.2.24 RC MinExposureX****Syntax**

```
RC ID←
```

Lower limit of the exposure multiplier.

**See Also**

[RC MaxExposureX](#) for upper limit of the exposure multiplier.

**Since**

Firmware 1.04.

**Examples**

```
>RC MinExposureX
OK:0:RC MinExposureX:1
```

---

#### 4.3.2.25 RC MaxExposureX

##### Syntax

```
RC MaxExposureX←
```

Upper limit of the exposure multiplier.

##### See Also

[RC MinExposureX](#) for lower limit of the exposure multiplier.

##### Since

Firmware 1.04.

##### Examples

```
>RC MaxExposureX
OK:0:RC MaxExposureX:50
```

#### 4.3.2.26 RC Speed

##### Syntax

```
RC Speed←
```

Lists the available speeds configured.

The response returns the number of speeds in the last parameter(4<sup>th</sup>) of the first line. Every line after that lists the speed ID and name.

##### Note

This command is only valid if the [RC InstrumentType](#) is 2 (Spectroradiometer).

##### Since

Firmware 1.17.

##### Examples

```
>RC Speed
OK:0:RC Speed:4
0,Slow
1,Normal
2,Fast
3,2x Fast
```

#### 4.3.2.27 RC MinSamplingRate

##### Syntax

```
RC MinSamplingRate←
```

Lower limit of the sampling rate reported in Hertz.

##### See Also

[RC MaxSamplingRate](#) for upper limit of sampling rate.

---

**Since**

Firmware 1.19.

**Examples**

```
>RC MinSamplingRate
OK:0:RC MinSamplingRate:200.0 Hz
```

**4.3.2.28 RC MaxSamplingRate****Syntax**

```
RC MaxSamplingRate←
```

Upper limit of the sampling rate reported in Hertz.

**See Also**

[RC MinSamplingRate](#) for lower limit of the sampling rate.

**Since**

Firmware 1.19.

**Examples**

```
>RC MaxSamplingRate
OK:0:RC MaxSamplingRate:1600.0 Hz
```

**4.3.3 Read Setup Commands**

The read setup command set is used to retrieve measurement properties from the instrument.

**Quick Reference**

Key	Description
Accessory	Current accessory
Filter	Current filters
Aperture	Current aperture
Mode	Current instrument mode
RangeMode	Current range mode
Range	Current range
ExposureMode	Current exposure mode
Exposure	Current exposure
SyncMode	Current sync mode
SyncFreq	Current sync frequency
ExposureX	Current exposure multiplier
MatrixMode	Current matrix mode <i>Deprecated as of Firmware 1.16.</i>
UserCalibMode	Current user calibration mode
Matrix	Current matrix calibration
Match	Current match calibration

Speed	Current speed
SamplingRate	Current sampling rate
MaxFreqFlickerSearch	Current Maximum Frequency flicker search
CMF	Current color matching function

#### 4.3.3.1 RS Accessory

##### Syntax

```
RS Accessory←
```

Retrieves the name of the current accessory.

##### Since

Firmware 1.04.

##### Examples

```
>RS Accessory
OK:0:RS Accessory:Standard
```

#### 4.3.3.2 RS Filter

##### Syntax

```
RS Filter←
```

Retrieves the names of the current filters as a comma separated value.

##### Since

Firmware 1.04.

##### Examples

```
>RS Filter
OK:0:RS Filter:ND-100-1,None,None
```

#### 4.3.3.3 RS Aperture

##### Syntax

```
RS Aperture←
```

Retrieves the name of the current aperture.

##### Since

Firmware 1.04.

##### Examples

```
>RS Aperture
OK:0:RS Aperture:5 deg
```

---

#### 4.3.3.4 RS Mode

##### Syntax

```
RS Mode←
```

Retrieves the name of the current instrument mode.

##### Since

Firmware 1.16.

##### Examples

```
>RS Mode  
OK:0:RS Mode:Colorimeter
```

#### 4.3.3.5 RS RangeMode

##### Syntax

```
RS RangeMode←
```

Retrieves the name of the current range mode.

##### Since

Firmware 1.04.

##### Examples

```
>RS RangeMode  
OK:0:RS RangeMode:Auto
```

#### 4.3.3.6 RS Range

##### Syntax

```
RS Range←
```

Retrieves the name of the current range. This is used only if the Range Mode retrieved by [RS RangeMode](#) is **Manual**.

##### Since

Firmware 1.04.

##### Examples

```
>RS Range  
OK:0:RS Range:A
```

#### 4.3.3.7 RS ExposureMode

##### Syntax

```
RS ExposureMode←
```

Retrieves the name of the current exposure mode.

---



**Since**

Firmware 1.04.

**Examples**

```
>RS ExposureMode
OK:0:RS ExposureMode:Auto
```

**4.3.3.8 RS Exposure****Syntax**

```
RS Exposure←
```

Retrieves the current exposure in msecs. This is used only if the Exposure Mode retrieved by [RS ExposureMode](#) is **Fixed**.

**Since**

Firmware 1.04.

**Examples**

```
>RS Exposure
OK:0:RS Exposure:1.000 msec
```

**4.3.3.9 RS SyncMode****Syntax**

```
RS SyncMode←
```

Retrieves the current sync mode.

**Since**

Firmware 1.04.

**Examples**

```
>RS SyncMode
OK:0:RS SyncMode:None
```

**4.3.3.10 RS SyncFreq****Syntax**

```
RS SyncFreq←
```

Retrieves the current sync frequency. This is used only if the Sync Mode retrieved by [RS SyncMode](#) is **Manual**.

**Since**

Firmware 1.04.

**Examples**

```
>RS SyncFreq
OK:0:RS SyncFreq:60.00 Hz
```

---

#### 4.3.3.11 RS ExposureX

##### Syntax

```
RS ExposureX←
```

Retrieves the current exposure multiplier.

##### Since

Firmware 1.04.

##### Examples

```
>RS ExposureX  
OK:0:RS ExposureX:1
```

#### 4.3.3.12 RS MatrixMode

##### Syntax

```
RS MatrixMode←
```

Retrieves the current matrix mode.

##### Since

Firmware 1.04.

**Deprecated** [RS MatrixMode](#) has been deprecated for [RS UserCalibMode](#) as of Firmware 1.16

##### Examples

```
>RS MatrixMode  
OK:0:RS Matrix:Disabled
```

#### 4.3.3.13 RS UserCalibMode

##### Syntax

```
RS UserCalibMode←
```

Reports the current user calibration mode set by the command [SM UserCalibMode](#).

##### Since

Firmware 1.16.

##### Examples

```
>RS UserCalibMode  
OK:0:RS UserCalibMode:None
```

#### 4.3.3.14 RS Matrix

##### Syntax

```
RS Matrix←
```

Retrieves the current calibration matrix ID. This is used only if the user calibration mode retrieved by [RS UserCalibMode](#) is **Matrix**.

---

**Since**

Firmware 1.04.

**Examples**

```
>RS Matrix  
OK:0:RS Matrix:0
```

**4.3.3.15 RS Match****Syntax**

```
RS Match←
```

Reports the current match calibration set ID which is selected by the command [SM Match](#). This is used only if the user calibration mode retrieved by [RS UserCalibMode](#) is **Match**.

**Since**

Firmware 1.16.

**Examples**

```
>RS Match  
OK:0:RS Match:0
```

**4.3.3.16 RS Speed****Syntax**

```
RS Speed←
```

Retrieves the name of the current speed.

**Note**

This command is only valid if the [RC InstrumentType](#) is 2 (Spectroradiometer).

**Since**

Firmware 1.17.

**Examples**

```
>RS Speed  
OK:0:RS Speed:Normal
```

**4.3.3.17 RS SamplingRate****Syntax**

```
RS SamplingRate←
```

Retrieves the current measurement sampling rate.

**Since**

Firmware 1.19.

**Examples**

```
>RS SamplingRate  
OK:0:RS SamplingRate:200.0 Hz
```

---

#### 4.3.3.18 RS MaxFreqFlickerSearch

##### Syntax

```
RS MaxFreqFlickerSearch←
```

Retrieves the current maximum frequency flicker search.

##### Since

Firmware 1.19.

##### Examples

```
>RS MaxFreqFlickerSearch  
OK:0:RS MaxFreqFlickerSearch:200.0 Hz
```

#### 4.3.3.19 RS CMF

##### Syntax

```
RS CMF←
```

Retrieves the current Color Matching Function Index.

##### Note

This command is only valid if the [RC InstrumentType](#) is 2 (Spectroradiometer).

##### Since

Firmware 1.26.

##### Examples

```
>RS CMF  
OK:0:RS CMF:0
```

note This

### 4.3.4 Measure Commands

The measure command set is used to capture data from the instrument.

##### Quick Reference

M Commands	Description
M	Capture a measurement
MT	Trigger a measurement
MA	Aborts a measurement in progress
MF	Measure the sync frequency

#### 4.3.4.1 M

##### Syntax

```
M←
```

Capture a measurement and returns after the measurement is completed or an error has occurred.

---

**Note**

This command does not require a command extension

**Since**

Firmware 1.04.

**See Also**

Refer 300 range [Response Codes](#) section.

**Examples**

```
>M
ER:-305:M:Light intensity too low or unmeasurable
>M
OK:0:M:No errors
```

**4.3.4.2 MT****Syntax**

```
MT←-
```

Trigger a measurement and returns immediately without waiting for the response.

**Note**

This is used for asynchronous capture of measurements.  
This command is reserved for future implementation.

**4.3.4.3 MA****Syntax**

```
MA←-
```

Aborts a measurement in progress.

**Note**

This is used for asynchronous capture of measurements.  
This command is reserved for future implementation.

**4.3.4.4 MF****Syntax**

```
MF←-
```

Measures the sync frequency of the light source.

**Note**

This command is reserved for future implementation.

## 4.4 Configure Commands

The configure command set is used to configure the instrument.

---

C Commands	Description
CC	<a href="#">Configure Calibration Commands</a>

#### 4.4.1 Configure Calibration Commands

The configure calibration command set is used to specify instrument calibration configuration properties.

##### Quick Reference

Key	Description
Matrix	Configures a new calibration matrix
Match	Configures a new match calibration

##### 4.4.1.1 CC Matrix

###### Syntax

```
CC Matrix AccIndex,MatrixID,MatrixName,R00,R01,R02,R10,R11,R12,R20,R21,R22←
```

Adds or updates a calibration matrix to the list of calibration matrices, where *AccIndex* is the index of the accessory, *MatrixID* is the index of the added or updated calibration matrix in the list, *MatrixName* is a description for the calibration matrix and *Rxy* are the elements of the matrix.

###### Since

Firmware 1.05.

###### Examples

```
>CC Matrix 0,0,Display Test,1.030e+00,-1.363e-02,-8.051e-03,-2.175e-02,1.-
072e+00,1.203e-02,5.340e-02,3.940e-03,1.058e+00
OK:0:CC Matrix:No errors
```

###### See Also

[RC MatrixCalib.](#)

##### 4.4.1.2 CC Match

###### Syntax

```
CC Match MatchID,MatchName,cfX,cfY,cfZ2←
```

Adds or updates a Match Calibration set to the list of match calibration sets, where *MatchID* is the index of the added or updated Match Calibration set in the list, *MatchName* is a description for the Match Calibration set and *cfX*, *cfY* and *cfZ* are the correction factor for the CIE tristimulus values X, Y, Z.

###### Since

Firmware 1.16.

###### Examples

```
>CC Match 0,Display Test,1.030e+00,-1.363e-02,-8.051e-03←
OK:0:CC Match:No errors
```

###### See Also

[RC MatchCalib.](#)

## 4.5 Response Codes

Following is a comprehensive list of response codes and their description.

### Note

Negative numbers are errors while positive numbers indicate warnings. 0 indicates no errors.

300 range is measurement errors.

500 range is general command errors

Error/Warning Code	Description
100	Light intensity too low for automatic sync
101	Cannot sync to constant light source
102	Cannot find sync, max limit selected
103	Sync level too low for reliable sync
-300	Invalid Sync mode
-301	Invalid Sync period
-302	Can not sync to light
-303	Light intensity is fluctuating
-304	Light intensity too low for range
-305	Light intensity too low or unmeasurable
-306	Light intensity too high for range
-307 to -330	Reserved
-331	Hardware malfunction
-332	Matrix version mismatch
-333	Invalid matrix index
-334	Uninitialized CIE tables
-335	Uninitialized CMF tables
-336	No Matrix exists for given ID
-500	Invalid command
-501	Reserved
-502	Reserved
-503	Reserved
-504	Reserved
-505	Duplicate Filter selection
-506	Index doesn't select an Accessory
-507	Index doesn't select a Filter
-508	Index not valid for Accessory
-509	Index not valid for Filter
-510	Index not valid for Filter
-511	Index not valid for Filter
-512	Invalid Range mode
-513	Invalid Range index
-514	Invalid Exposure Multiplier
-515	Index doesn't select an Aperture
-516	Reserved

-517	Reserved
-518	Invalid Exposure Mode
-519	Invalid Exposure value
-520	Reserved
-521	Invalid Sync Mode
-522	Invalid User Sync Frequency
-523	Reserved
-524	Reserved
-525	Reserved
-526	Reserved
-527	Reserved
-528	Reserved
-529	Reserved
-530	Reserved
-531	Reserved
-532	Reserved
-533	Reserved
-534	Reserved
-535	Reserved
-536	Reserved
-537	Reserved
-538	Reserved
-539	Reserved
-540	Reserved
-541	Reserved
-542	Reserved
-543	Reserved
-544	Reserved
-545	Reserved
-546	Reserved
-547	Reserved
-548	Reserved
-549	Reserved
-550	Reserved
-551	Reserved
-552	Invalid Matrix Mode
-553	Invalid Matrix ID
-555	Invalid Matrix name/description
-556	Error saving Matrix to FLASH
-557	Invalid Match ID
-558	Invalid Match name/description
-559	Error saving Match to FLASH
-560	Invalid User Calibration Mode





## Chapter 5

# Deprecated List

### Page [Remote Commands](#)

[SM MatrixMode](#) has been deprecated for [SM UserCalibMode](#) as of Firmware 1.16

[RM MatrixMode](#) has been deprecated for [RM UserCalibMode](#) as of Firmware 1.16

[RC MatrixMode](#) has been deprecated for [RC UserCalibMode](#) as of Firmware 1.16

[RC MatrixCalibration](#) has been deprecated for [RC MatrixCalib](#) as of Firmware 1.16

[RS MatrixMode](#) has been deprecated for [RS UserCalibMode](#) as of Firmware 1.16